

RNA binding protein site analysis with CLIP

This document describes a protocol for preparing performing an analysis of RNA binding protein (RBP) interaction sites on RNA. This protocol is based on the single-end, monitored eCLIP, and irCLIP2 protocols described previously (PMID: 28766298, 30252095, 31610236). Compared to these protocols, I use modifications at various steps during processing, as well as in the sequencing adapter configuration. Libraries made with this protocol can be sequenced on Illumina platforms, such as NextSeq, using single- or paired-end reads. The UMI strategy used is based on PMID: 30001700.

Quick links

- 1. Cell culture and harvest
 - 1.1 Reagents and materials
 - 1.2 Solution recipes
 - 1.3 Protocol
- 2. Cell lysis, immunoprecipitation, and ligation
 - 2.1 Reagents and materials
 - 2.2 Solution recipes
 - 2.3 Oligos
 - 2.4 Protocol
- 3. SDS-PAGE and RNA isolation
 - 3.1 Reagents and materials
 - 3.2 Solution recipes
 - 3.3 Protocol
- 4. Phosphorylation, biotin capture, reverse transcription, and linker ligation
 - 4.1 Reagents and materials
 - 4.2 Solution recipes
 - 4.3 Oligos
 - 4.4 Protocol
- 5. Library amplification and size selection
 - 5.1 Reagents and materials
 - 5.2 Solution recipes
 - 5.3 Oligos
 - 5.4 Protocol
- Notes
- Reagent list

1. Cell culture and harvest

Use whatever cell culture methods you are comfortable with here. This stage of the protocol is only given as an example of a process that is validated to work and is by no means the only way to achieve a successful result in this section.

1.1 Reagents and materials

- 15mL tubes (VWR, CAT#89093-186)
- 50mL tubes (VWR, CAT#89093-190)
- 5, 10, 25mL pipettes (multiple vendors)
- Centrifuge with 15mL tube buckets (multiple vendors)
- TrypLE express (Thermo Scientific, CAT#12605028)
- DPBS, no calcium or magnesium (Thermo Scientific, CAT#14190144)
- Stratalinker UV Crosslinker, 254nm (Stratagene)
- Dry ice

1.2 Solution recipes

- Cell culture medium (provided only as an example)
 - DMEM + GlutaMAX (500mL)
 - FBS (50mL)

1.3 Protocol

This protocol assumes you have your cells growing in 15cm dishes. For other dish sizes, simply scale the reagents for your purpose. Because you need to crosslink, I recommend only processing one to two plates at a time.

1. Label a 15mL tube for each sample you plan to harvest.
2. Label a 50mL tube with 'TE', another with 'Media', and a last with 'PBS'.
3. To the 'TE' 50mL tube, add 8mL of TrypLE express per sample you plan to harvest.
4. To the 'PBS' 50mL tube, add 5mL of DPBS per sample you plan to harvest.
5. Remove the cell plate from the incubator and pipette the media the cells are growing in into the 50mL tube labeled 'Media'.
6. Carefully add 10mL of room temperature (+21C) DPBS to the plate, swirl, and discard the DPBS rinse to the waste.
7. Carefully add 5mL of room temperature DPBS to the plate and place on an ice block with some water on top in the Stratalinker (make sure the plate is flat). Crosslink at 200mJ/cm² (see **Note 1**).
8. Remove the DPBS and add 8mL of TrypLE express to the plate and place back into the incubator.
9. Once the cells have started to detach, remove the plate from the incubator.
10. With 4mL of media from the 'Media' tube, rinse and collect the detached cells and transfer to the appropriate 15mL tube.
11. Centrifuge the tube with the cells for 3-minutes at 250g.
12. Dump the supernatant to waste.
13. Add 5mL of DPBS from the 'PBS' tube to the tube with the cell pellet. There is no need to pipette mix here.
14. Centrifuge the tube for 30-seconds at 250g.

15. Dump the supernatant to waste and spirate the remaining liquid to leave as little as possible.
16. Place the cell pellet on dry ice until frozen and keep at -80C until use.

2. Cell lysis, immunoprecipitation, and ligation

You should use an antibody that you have validated has good selectivity for your RBP of interest.

2.1 Reagents and materials

- 1.5mL or 2.0mL Safe-Lock tubes (Thermo Scientific, CAT#05-402-25 or CAT#05-402-7)
- Benchtop centrifuge with holder for microcentrifuge tubes (multiple vendors)
- ThermoMixer with a 2mL block (Eppendorf)
- Magnetic rack for 1.5mL or 2.0mL tubes (Promega, CAT#Z5342)
- RNaseZAP wipes (Thermo Scientific, CAT#AM9788)
- TURBO DNase, 2U/uL (Thermo Scientific, CAT#AM2238)
- RNaseIn, 40U/uL (Promega, CAT#N2511)
- 1M TrisCl pH 7.5 (Thermo Scientific, CAT#15567027)
- 5M NaCl (Thermo Scientific, CAT#AM9760G)
- 0.5M EDTA (Thermo Scientific, CAT#15575020)
- NP-40 (Sigma, CAT#I8896) - NP-40 seems to no longer be available, so we use IGEPAL CA-630 as a substitute
- Deoxycholate (DOC) (Sigma, CAT#D6750)
- 1M MgCl2 (Thermo Scientific, CAT#AM9530G)
- cOmplete protease inhibitor tablets, EDTA free (Sigma, CAT#11836170001)
- Qubit assay tubes (Thermo Scientific, CAT#Q32856)
- Qubit RNA HS Assay Kit (Thermo Scientific, CAT#Q32852)
- RNase 1, 100U/uL (Thermo Scientific, CAT#AM2294)
- RNA clean water (Thermo Scientific, CAT#10977023)
- FastPrep 5G instrument (MP Biomedicals)
- Lysing matrix Y 2mL tubes (MP Biomedicals, CAT#116960050)
- NuPAGE LDS Sample Buffer (4X) (Thermo Scientific, CAT#NP0007)
- Sucrose (Thermo Scientific, CAT#15503022)
- 20% SDS (Thermo Scientific, CAT#BP1311)
- Protein A/G magnetic beads (Thermo Scientific, CAT#88802)
- Bovine serum albumin, fraction V (Sigma, CAT#A2153)
- Dithiothreitol (Bio-Rad, CAT#1610611)
- PEG-8000 (Promega, CAT#V3011)
- QuickCIP (NEB, CAT#M0525L)
- T4 PNK (NEB, CAT#M0201L)
- T4 RNA Ligase 2, truncated KQ (NEB, CAT#M0373L)
- 10mM ATP (NEB, CAT#P0756S)
- T4 RNA ligase 1 (NEB, CAT#M0437M)

2.2 Solution recipes

- 10% (v/v) NP-40 - 1mL in 9mL of water
- 10% (w/v) DOC - 1g in 10mL of water
- 50% (w/v) sucrose - 5g in 10mL of water (filter and keep at +4C)
- 10X cOmplete protease inhibitor stock - 1 tablet in 1mL of water
- 50% (w/v) PEG-8000 - 12.5g in 25mL of water
- 0.2M of dithiothreitol (DTT) - 15mg in 500uL of water, prepare fresh and keep on ice
- 1M of dithiothreitol (DTT) - 75mg in 500uL of water, prepare fresh and keep on ice
- Diluted RNase 1 - 5uL of 100U/uL stock in 1mL of water (now 0.5U/uL)
- Lysis buffer (need 1mL per sample)
 - 50mM Tris-Cl pH 7.5 (50uL of 1M stock)
 - 150mM NaCl (30uL of 5M stock)
 - 10% (w/v) sucrose (200uL of 50% stock)
 - 0.5% (v/v) NP-40 (50uL of 10% stock)
 - 0.5% (v/v) DOC (50uL of 10% stock)
 - 0.1% (v/v) SDS (5uL of 20% stock)
 - 0.5X cOmplete protease inhibitor (50uL of 10X stock)
 - 20U/mL TURBO DNase (10uL of 2U/uL stock)
 - 200U/mL SUPERaseIn (5uL of 40U/uL stock)
 - water to 1mL (550uL of pre-chilled water)
- Rinse buffer A (make 50mL and keep at +4C)
 - 50mM Tris-Cl pH 7.5 (2.5mL of 1M stock)
 - 100mM NaCl (1mL of 5M stock)
 - 0.5% (v/v) NP-40 (2.5mL of 10% stock)
 - 0.5% (v/v) DOC (2.5mL of 10% stock)
 - 0.1% (v/v) SDS (250uL of 20% stock)
 - water to 50mL (42mL of pre-chilled water)
- Rinse buffer B (make 50mL and keep at +4C)
 - 50mM Tris-Cl pH 7.5 (2.5mL of 1M stock)
 - 100mM NaCl (1mL of 5M stock)
 - 0.5% (v/v) NP-40 (2.5mL of 10% stock)
 - 0.5% (v/v) DOC (2.5mL of 10% stock)
 - 0.1% (v/v) SDS (250uL of 20% stock)
 - 0.1mg/mL BSA (5mg)
 - water to 50mL (42mL of pre-chilled water)
- Rinse buffer C (make 50mL and keep at +4C)
 - 50mM Tris-Cl pH 7.5 (2.5mL of 1M stock)
 - 1M NaCl (10mL of 5M stock)
 - 5mM EDTA (500uL of 0.5M stock)
 - 0.5% (v/v) NP-40 (2.5mL of 10% stock)
 - 0.5% (v/v) DOC (2.5mL of 10% stock)
 - 0.1% (v/v) SDS (250uL of 20% stock)
 - water to 50mL (32mL of pre-chilled water)

- Rinse buffer D (make 50mL and keep at +4C)
 - 50mM Tris-Cl pH 7.5 (2.5mL of 1M stock)
 - 10mM MgCl₂ (500uL of 1M stock)
 - water to 50mL (47mL of pre-chilled water)
- PNK buffer 1 (make 1mL per sample)
 - 50mM Tris-Cl pH 7.5 (50uL of 1M stock)
 - 10mM MgCl₂ (10uL of 1M stock)
 - 1mM DTT (5uL of 200mM stock)
 - water to 1mL (935uL)
- RNL buffer 1 (make 1mL per sample)
 - 50mM Tris-Cl pH 7.5 (50uL of 1M stock)
 - 10mM MgCl₂ (10uL of 1M stock)
 - 1mM DTT (5uL of 200mM stock)
 - 20% (v/v) PEG-8000 (400uL of 50% stock)
 - water to 1mL (535uL)
- PNK buffer 2 (make 1mL per sample)
 - 50mM Tris-Cl pH 7.5 (50uL of 1M stock)
 - 10mM MgCl₂ (10uL of 1M stock)
 - 1mM DTT (5uL of 200mM stock)
 - 1mM ATP (100uL of 10mM stock)
 - water to 1mL (835uL)
- RNL buffer 2 (make 1mL per sample)
 - 50mM Tris-Cl pH 7.5 (50uL of 1M stock)
 - 10mM MgCl₂ (10uL of 1M stock)
 - 1mM DTT (5uL of 200mM stock)
 - 1mM ATP (100uL of 10mM stock)
 - 20% (v/v) PEG-8000 (400uL of 50% stock)
 - water to 1mL (435uL)

2.3 Oligos

- Ribo_3DNALinker_v4 (100nmol scale, RNase-free HPLC purified)
 - /5rApp/NNNGTCNNNTAGNNNAGATCGGAAGAGCACACGTCT/3ddC/
 - reconstitute to 20uM and aliquot. Store at -20C.
- Clip_5RNALinker_v4 (100nmol scale, RNase-free HPLC purified)
 - /5BiotinTEG/rCrArGrUrCrCrGrArCrGrArUrC
 - reconstitute to 20uM and aliquot. Store at -20C.

2.4 Protocol

Read this protocol before starting to make sure you are prepared for the steps as they come. All of the buffers should be cold (+4C). Pre-chill a benchtop centrifuge to +4C before starting. Pre-heat a Thermomixer to +37C before starting.

1. Prepare the antibody-conjugated beads (all rinses are 5-minutes at room temperature (+21C)):

1. Wash 125uL of protein A/G beads with 500uL of Rinse buffer A in a 2mL tube. Discard the supernatant.
2. Repeat the previous step two more times for a total of 3 rinses.
3. Reconstitute the washed beads in 125uL of Rinse buffer A and transfer 25uL of this to a fresh 2mL tube (set this tube aside on ice).
4. Add 10ug of your target antibody to the remaining 100uL of beads in Rinse buffer A. Rotate slowly at room temperature for 1-hour (see **Note 2**).
5. Rinse the conjugated beads twice with 500uL of Rinse buffer B (with BSA) (see **Note 3**).
6. Rinse the conjugated beads three times with Rinse buffer A.
7. Reconstitute the beads in 100uL of Rinse buffer A until use. Keep on ice.
2. Prepare your lysate and perform IP:
 1. Slowly thaw cell pellet at room temperature until a small amount of frozen material remains.
 2. Add 500uL of Lysis buffer to the pellet and pipette mix. Incubate on ice for 5-minutes.
 3. Transfer the lysate to a Lysing matrix Y tube and disrupt on a FastPrep (6M/s, 45 seconds, 1 cycle) (see **Note 4**).
 4. Centrifuge the lysate at 12,000g for 5-minutes and recover the supernatant to a fresh 2mL tube (see **Note 5**).
 5. Take 5uL of the lysate and dilute to 50uL with water in a fresh 1.5mL tube.
 6. Measure the RNA concentration of 5uL of this diluted sample using the Qubit RNA assay. Calculate the RNA concentration in the original sample.
 7. Transfer a volume equivalent to 20ug of RNA to a fresh 2mL tube and bring to 400uL total volume with Lysis buffer.
 8. Add 1 Unit of RNase 1 per 10ug of measured RNA in the sample (4uL of Diluted RNase 1 for 20ug sample). Pipette mix and incubate at +37C for 5-minutes (see **Note 6**).
 9. Transfer the supernatant to the tube containing 25uL of rinsed protein A/G beads prepared above. Incubate with slow rotation at +4C for 30-minutes.
 10. Recover the supernatant lysate and combine with the 100uL of conjugated beads prepared above. Incubate with slow rotation at +4C for 2-hours (see **Note 7**).
 11. Recover 10uL (2%) of the lysate-bead mix to a fresh 2mL tube (**SampleInput**) and another 10uL to a different 2mL tube (**ImagingInput**). Bring the volume of each to 20uL with Lysis buffer.
 12. Preheat the ThermoMixer to +70C.
 13. To the Input samples, add 8uL of NuPAGE LDS Sample Buffer (4X) and 4uL of 1M DTT and pipette mix.
 14. Incubate the Input samples at +70C with mixing at 1,000rpm for

- 10-minutes.
15. Recover the supernatant from the eluted Input samples to fresh 1.5mL tubes and store at -20C until SDS-PAGE.
 16. Reduce ThermoMixer heat to +37C.
 17. Separate the remaining supernatant from the Ip beads and transfer to a fresh 2mL tube. This supernatant can be saved as a control to check performance of the Ip.
 18. Rinse the remaining beads twice with 500uL of Rinse buffer C. Discard the supernatant each time.
 19. Rinse the beads twice in 500uL of Rinse buffer A. Discard the supernatant each time.
 20. Rinse the beads twice with 500uL of Rinse buffer D. Discard the supernatant each time.
 21. Rinse the beads once with 500uL of PNK buffer 1. Discard the supernatant.
3. Perform on-bead dephosphorylation and ligation:
1. Reconstitute the rinsed beads in 50uL of PNK buffer 1. Add 2uL of QuickCIP and 2uL of T4 PNK and pipette mix. Incubate at +37C for 30-minutes with mixing at 1,000rpm.
 2. Remove the supernatant and rinse once with 500uL of Rinse buffer C. Discard the supernatant.
 3. Rinse the beads twice in 500uL of Rinse buffer A. Discard the supernatant each time.
 4. Rinse the beads twice with 500uL of Rinse buffer D. Discard the supernatant each time.
 5. Rinse the beads once in 500uL of RNL buffer. Discard the supernatant.
 6. Reconstitute the rinsed beads in 45uL of RNL buffer 1. Add 5uL of Ribo_3DNALinker_v4 followed by 1uL of T4 RNA Ligase 2 truncated KQ (200 Units). Pipette mix and incubate at room temperature for 4-hours with gentle mixing at 1,000rpm.
 7. Remove the supernatant and rinse once with 500uL of Rinse buffer C. Discard the supernatant.
 8. Rinse the beads twice in 500uL of Rinse buffer A. Discard the supernatant each time.
 9. Rinse the beads twice with 500uL of Rinse buffer D. Discard the supernatant each time.
 10. Rinse the beads once in 500uL of PNK buffer 2. Discard the supernatant.
 11. Reconstitute the rinsed beads in 50uL of PNK buffer 2. Add 2uL of T4 PNK (20 Units) and pipette mix. Incubate at +37C for 30-minutes with mixing at 1,000rpm.
 12. Rinse the beads twice in 500uL of Rinse buffer A. Discard the supernatant each time.
 13. Rinse the beads twice with 500uL of Rinse buffer D. Discard the supernatant each time.
 14. Rinse the beads once in 500uL of RNL buffer 2. Discard the super-

- natant.
15. Reconstitute the rinsed beads in 45uL of RNL buffer 2. Add 5uL of Clip_5RNALinker_v4 and 2uL of T4 RNA ligase 1 (60 Units). Pipette mix and incubate at room temperature overnight with gentle mixing at 1,000rpm.
 16. Preheat the ThermoMixer to +70C.
 17. Remove the supernatant and rinse the beads twice in 500uL of Rinse buffer A. Discard the supernatant each time.
 18. Rinse the beads twice with 500uL of Rinse buffer D. Discard the supernatant each time.
 19. Reconstitute the beads in 25uL of Lysis buffer (**SampleIp**) and transfer 5uL (20%) of the mix to a fresh 2mL tube containing 15uL of Lysis buffer (**ImagingIp**).
 20. Add 8uL of NuPAGE LDS Sample Buffer (4X) and 4uL of 1M DTT to each tube and pipette mix. Incubate the Ip samples at +70C with mixing at 1,000rpm for 10-minutes.
 21. Recover the supernatant from the eluted Ip samples to fresh 1.5mL tubes and store at -20C until SDS-PAGE.

3. SDS-PAGE and RNA isolation

3.1 Reagents and materials

- 1.5mL or 2.0mL Safe-Lock tubes (Thermo Scientific, CAT#05-402-25 or CAT#05-402-7)
- Benchtop centrifuge with holder for microcentrifuge tubes (multiple vendors)
- ThermoMixer with a 2mL block (Eppendorf)
- RNaseZAP wipes (Thermo Scientific, CAT#AM9788)
- PageRuler Prestained Protein Ladder, 10 to 180 kDa (Thermo Scientific, CAT#26616)
- NuPAGE 4 to 12%, Bis-Tris, 1.5 mm, Mini Protein Gel, 10-well (Thermo Scientific, CAT#NP0335BOX)
- NuPAGE MOPS SDS Running Buffer (20X) (Thermo Scientific, CAT#NP0001)
- XCell SureLock Mini-Cell Electrophoresis System (Thermo Scientific)
- 1M TrisCl pH 7.5 (Thermo Scientific, CAT#15567027)
- 1M MgCl₂ (Thermo Scientific, CAT#AM9530G)
- RNA clean water (Thermo Scientific, CAT#10977023)
- Dithiothreitol (Bio-Rad, CAT#1610611)
- NuPAGE Transfer Buffer (20X) (Thermo Scientific, CAT#NP0006)
- Amersham Protran Premium Western blotting membranes, nitrocellulose (Sigma, CAT#GE10600008)
- PVDF Transfer Membrane, 0.45um (Thermo Scientific, CAT#88518)
- Methanol, HPLC-grade (Sigma, CAT#34860-4L-R)
- TrueBlot Anti-Rabbit IgG HRP (Cedarlane, CAT#18-8816-31)

- DPBS, no calcium, no magnesium (Thermo Scientific, CAT#14190144)
- Powdered milk (grocery store)
- Tween 20 (Sigma, CAT#P9416)
- Tris base (Sigma, CAT#93362)
- NaCl (Sigma, CAT#71376)
- 12N HCl (Sigma, CAT#H1758)
- Pierce ECL Western Blotting Substrate (Thermo Scientific, CAT#32109)
- CL-XPosure Film (Thermo Scientific, CAT#34089)
- 0.5M EDTA (Thermo Scientific, CAT#15575020)
- 5M NaCl (Thermo Scientific, CAT#AM9760G)
- Proteinase K Solution (20 mg/mL) (Thermo Scientific, CAT#AM2546)
- Urea (Sigma, CAT#U5378)
- Razorblades (any will do)
- Acid-Phenol:Chloroform, pH 4.5 (with IAA, 125:24:1) (Thermo Scientific, CAT#AM9720)
- Phenol:Chloroform:Isoamyl Alcohol 25:24:1, Saturated with 10mM Tris, pH 8.0, 1mM EDTA (Sigma, CAT#P3803-100ML)
- Phaselock gel (heavy) (VWR, CAT#CA10847-802)
- RNA Clean & Concentrator-5 kit (Cedarlane, CAT#R1015)
- Absolute ethanol (wherever your institute normally procures this from)

3.2 Solution recipes

- Diluted marker - 2uL PageRuler Prestained Protein Ladder, 3uL NuPAGE LDS Sample Buffer (4X), 6uL water
- Rinse buffer D (make 50mL)
 - 50mM Tris-Cl pH 7.5 (2.5mL of 1M stock)
 - 10mM MgCl₂ (500uL of 1M stock)
 - water to 50mL (47mL of pre-chilled water)
- 10X TBS (make 1L)
 - 200mM Tris base (24g)
 - 1.5M NaCl (88g)
 - water to 500mL
 - adjust to pH 7.6 with 12N HCl
 - water to 1L
- 1X TBST (make 1L)
 - 1X TBS (100mL)
 - water to 1L (900mL)
 - Tween 20 (1mL)
- 5% milk in TBST (make 100mL)
 - powdered milk (5g)
 - TBST to 100mL
- Proteinase K buffer 1 (make 10mL)
 - 100mM Tris-Cl pH 7.5 (1mL of 1M stock)
 - 10mM EDTA (200uL of 0.5M stock)
 - 50mM NaCl (100uL of 5M stock)

- water to 10mL (8.7mL)
- Proteinase K buffer 2 (make 10mL)
 - 100mM Tris-Cl pH 7.5 (1mL of 1M stock)
 - 10mM EDTA (200uL of 0.5M stock)
 - 50mM NaCl (100uL of 5M stock)
 - 7M urea (4.2g)
 - water to 10mL

3.3 Protocol

Read this protocol before starting to make sure you are prepared for the steps as they come. All of the buffers should be cold (+4C).

1. For SDS-PAGE (you will have 2 gels, one for Samples, one for Imaging):
 1. Prepare the gels and pre-run for 10-minutes at 150V.
 2. Load samples on the gel with 10uL of the PageRuler marker on each end, alternating samples with Diluted marker (10uL each lane).
 3. For all samples, load the entire 30uL. Run at 150V for 75-minutes or until the dye front reaches the bottom.
2. Perform the transfer:
 1. For the imaging gel, activate a PVDF membrane in methanol and transfer to 1X Transfer buffer.
 2. For the sample gel, transfer nitrocellulose membrane to 1X Transfer buffer.
 3. Assemble transfer stacks, from bottom to top (black side (negative) of stack holder on bottom): (negative) 1x sponge – 2x filter paper – gel – membrane – 2x filter paper – 1x sponge (positive).
 4. Transfer overnight at 30V.
 5. The next day, remove the sample gel membrane, rinse once with DPBS, then wrap in Saran wrap and store at -20C.
3. To develop the imaging gel membrane:
 1. Block membrane in 5% milk in TBST at room temperature for 30-minutes.
 2. Probe with primary antibody at appropriate concentration (typically 0.2ug/mL) in 5% milk in TBST, at room temperature for 1-hour.
 3. Wash three times with TBST, for 5-minutes each.
 4. Probe with secondary antibody: 1:4000 Rabbit TrueBlot HRP in 5% milk in TBST and incubate at room temperature for 1-3 hours.
 5. Wash three times with TBST, for 5-minutes each.
 6. Develop with ECL, 30-seconds to 5-minutes, and image with western blot film.
4. Cutting the sample gel membrane:
 1. Place membrane on clean glass/metal surface.
 2. Using a fresh razor blade, cut desired lane from the RBP band to a region ~75kDa above your target.
 3. Slice the membrane piece into ~1–2 mm pieces.

4. Transfer slices to a 2mL tube on ice. Remember to also cut the same region from your SampleInput lane.
5. RNA extraction:
 1. To the membrane pieces, add 180uL of Proteinase K buffer 1 and 20uL of Proteinase K (20mg/mL). Incubate at +37C for 20-minutes in a ThermoMixer at 1,000rpm.
 2. Add 180uL of Proteinase K buffer 2 and 20uL of Proteinase K (20mg/mL) and incubate for an additional 20-minutes at +37C with mixing at 1,000rpm.
 3. Add 400uL of Acid phenol chloroform to the SampleInput tubes and 400uL of Phenol chloroform (non-acidic) to the SampleIp tubes and mix well. Incubate at +37C for 5-minutes.
 4. Transfer liquid (no membrane pieces) to a Phaselock gel (heavy) tube and incubate in ThermoMixer at 1,000 rpm, +37C for 5-minutes.
 5. Centrifuge at 13,000g for 15-minutes at room temperature.
 6. Transfer the aqueous (top) layer to a new 15 mL (or at least 3 mL volume) conical tube.
 7. Add 2 volumes RNA binding buffer from the RNA Clean and Concentrator kit (typically ~800uL).
 8. Add equal volume 100% ethanol and mix (typically ~1200uL).
 9. Transfer 750uL of mixed sample to a Zymo-Spin column.
 10. Centrifuge for 30-seconds at 12,000g and discard the flow-through.
 11. Repeat spins by reloading an additional 750uL volume until all sample has been spun through the column.
 12. Add 400uL of RNA Prep Buffer, centrifuge for 30-seconds at 12,000g, and discard the flow-through.
 13. Add 700uL of RNA Wash Buffer, centrifuge for 30-seconds at 12,000g, and discard the flow-through.
 14. Add 400uL of RNA Wash Buffer, centrifuge for 30-seconds at 12,000g, and discard the flow-through.
 15. Centrifuge for additional 2-minutes at 12,000g.
 16. Transfer column to a new 1.5mL tube, add 25uL of water to the column, incubate for 1-minute, and then centrifuge for 30-seconds at 12,000g.
 17. Samples can be stored at -80C.

4. Phosphorylation, biotin capture, reverse transcription, and linker ligation

4.1 Reagents and materials

- 1.5mL or 2.0mL Safe-Lock tubes (Thermo Scientific, CAT#05-402-25 or CAT#05-402-7)
- PCR microtubes (any vendor is fine)
- Benchtop centrifuge with holder for microcentrifuge tubes (multiple vendors)

- ThermoMixer with a 2mL block (Eppendorf)
- Magnetic rack for 1.5mL or 2.0mL tubes (Promega, CAT#Z5342)
- PCR Thermocycler (multiple vendors)
- RNaseZAP wipes (Thermo Scientific, CAT#AM9788)
- RNA clean water (Thermo Scientific, CAT#10977023)
- Dithiothreitol (Bio-Rad, CAT#1610611)
- PEG-8000 (Promega, CAT#V3011)
- QuickCIP (NEB, CAT#M0525L)
- T4 PNK (NEB, CAT#M0201L)
- T4 RNA Ligase 2, truncated KQ (NEB, CAT#M0373L)
- 10mM ATP (NEB, CAT#P0756S)
- Dynabeads MyOne Silane (Thermo Scientific, CAT#37002D)
- RLT buffer (Qiagen, CAT#79216)
- TURBO DNase, 2U/uL (Thermo Scientific, CAT#AM2238)
- RNaseIn, 40U/uL (Promega, CAT#N2511)
- 1M TrisCl pH 7.5 (Thermo Scientific, CAT#15567027)
- 5M NaCl (Thermo Scientific, CAT#AM9760G)
- 0.5M EDTA (Thermo Scientific, CAT#15575020)
- 1M MgCl₂ (Thermo Scientific, CAT#AM9530G)
- Dynabeads MyOne Streptavidin C1 (Thermo Scientific, CAT#65001)
- Absolute ethanol (wherever your institute purchases this from)
- NP-40 (Sigma, CAT#I8896) - NP-40 seems to no longer be available, so we use IGEPAL CA-630 as a substitute
- Exonuclease I, 20U/uL (E. coli) (NEB, CAT#M0293S) - comes with 10X Exo1 buffer
- 1N NaOH (Sigma, CAT#S2770-100ML)
- 1N HCl (Sigma, CAT#H9892-100ML)
- T4 RNA ligase 1 (NEB, CAT#M0437M)
- SuperScript IV Reverse Transcriptase (Thermo Scientific, CAT#18090050)
- this includes 5X SSIV buffer and 0.1M DTT
- Deoxynucleotide (dNTP) Solution Mix, 10mM (NEB, CAT#N0447S)

4.2 Solution recipes

- 50% (w/v) PEG-8000 - 12.5g in 25mL of water
- 0.2M of dithiothreitol (DTT) - 15mg in 500uL of water, prepare fresh and keep on ice
- 80% ethanol (10mL of ethanol in 40mL of water)
- 10% (v/v) NP-40 - 1mL in 9mL of water
- 2X PNK buffer 1 (make 1mL)
 - 100mM Tris-Cl pH 7.5 (100uL of 1M stock)
 - 20mM MgCl₂ (20uL of 1M stock)
 - 2mM DTT (10uL of 200mM stock)
 - water to 1mL (870uL)
- 2X RNL buffer 1 (make 1mL)
 - 100mM Tris-Cl pH 7.5 (100uL of 1M stock)

- 20mM MgCl₂ (20uL of 1M stock)
- 2mM DTT (10uL of 200mM stock)
- 40% (w/v) PEG-8000 (800uL of 50% stock)
- water to 1mL (70uL)
- Streptavidin solution A (make 10mL)
 - 100mM NaOH (1mL from 1N stock)
 - 50mM NaCl (100uL from 5M stock)
 - 0.1% NP-40 (100uL from 10% stock)
- Streptavidin solution B (make 10mL)
 - 100mM NaCl (200uL of 5M stock)
 - 0.1% NP-40 (100uL from 10% stock)
- 2X B&W buffer (make 10mL)
 - 10mM TrisCl pH 7.5 (100uL of 1M stock)
 - 1mM EDTA (20uL of 0.5M stock)
 - 2M NaCl (4mL of 5M stock)
 - 0.2% NP-40 (200uL of 10% stock)
- 1X B&W buffer (make 5mL) - 2.5mL of 2X B&W buffer with 2.5mL of water
- 10uM Ribo_RTPrimer_v4 - 10uL of 100uM stock in 90uL of water
- 2X RT master mix (make 10uL per sample)
 - 2X SSIV buffer (4uL of 5X stock)
 - 1mM dNTPs (1uL of 10mM stock)
 - 10mM DTT (1uL of 0.1M stock)
 - 4U/uL RNaseIn (1uL of 40U/uL stock)
 - 20U/uL SSIV (1uL of 200U/uL stock)
 - water (2uL)
- RNL buffer 2 (make 1mL)
 - 50mM Tris-Cl pH 7.5 (50uL of 1M stock)
 - 10mM MgCl₂ (10uL of 0.5M stock)
 - 1mM DTT (5uL of 200mM stock)
 - 1mM ATP (100uL of 10mM stock)
 - 20% (v/v) PEG-8000 (400uL of 50% stock)
 - water to 1mL (435uL)
- 3' linker mix - mix equal parts of 40uM stocks of Clip_3DNALinker2_v4a and Clip_3DNALinker2_v4b

4.3 Oligos

- Ribo_3DNALinker_v4 (100nmol scale, RNase-free HPLC purified)
 - /5rApp/NNNGTCNNNTAGNNNAGATCGGAAGAGCACACGTCT/3ddC/
 - reconstitute to 20uM and aliquot. Store at -20C.
- Ribo_RTPrimer_v4 (100nmol scale, desalting purification)
 - AGACGTGTGCTCTTCCGATCT
 - reconstitute at 100uM. Store at -20C.
- Clip_3DNALinker2_v4a (100nmol scale, HPLC purified)
 - /5Phos/NNNGTANNNTCGNNNAGATCGGAAGAGCGTCGTG/3ddC/

- reconstitute to 40uM and aliquot. Store at -20C.
- Clip_3DNALinker2_v4b (100nmol scale, HPLC purified)
 - /5Phos/NNNACTNNNGATNNNAGATCGGAAGAGCGTCGTG/3ddC/
 - reconstitute to 40uM and aliquot. Store at -20C.

4.4 Protocol

The starting point for this stage of the protocol is the 10uL purified RNA sample from above.

1. First we need to 3' dephosphorylate the SampleInput RNA:
 1. Add 25uL of 2X PNK buffer 1 to the input sample only and then add 2uL of QuickCIP and 1uL of T4 PNK (10 Units of each) and pipette mix. Incubate for 30-minutes at +37C.
 2. In a separate tube, rinse 20uL of MyONE Silane beads (per sample) with 500uL of RLT buffer. Discard the supernatant.
 3. Resuspend the beads in 365uL of RLT buffer and add to the input sample. Pipette mix.
 4. Add 10uL of 5M NaCl to the sample and pipette mix.
 5. Add 400uL of ethanol and rotate the tubes to mix for 15-minutes at room temperature (+21C).
 6. Place tubes on magnetic rack and discard the supernatant.
 7. Rinse the beads three times with 800uL of 80% ethanol. Discard supernatants to the waste.
 8. Add 40uL of water to the beads and incubate at room temperature for 5-minutes.
 9. Recover 20uL of this sample and store at -20C as a backup. Keep the other 20uL with the beads.
2. Now perform 3' adapter ligation of the SampleInput RNA:
 1. Add 25uL of 2X RNL buffer 1 to the bead sample, and then add 5uL of Ribo_3DNALinker_v4 followed by 1uL of T4 RNA Ligase 2 truncated KQ (200 Units). Pipette mix and incubate at room temperature for 4-hours with gentle rotation.
 2. Resuspend the beads in 365uL of RLT buffer and add to the input sample. Pipette mix.
 3. Add 10uL of 5M NaCl to the sample and pipette mix.
 4. Add 400uL of ethanol and rotate the tubes to mix for 15-minutes at room temperature (+21C).
 5. Place tubes on magnetic rack and discard the supernatant.
 6. Rinse the beads three times with 800uL of 80% ethanol. Discard supernatants to the waste.
 7. Add 10uL of water to the beads and incubate at room temperature for 5-minutes.
 8. Recover 8uL of the sample from the beads and keep at on ice or at -20C until reverse transcription.
3. SampleIp RNA can now be enriched on biotin beads to remove unligated

fragments:

1. Wash 10uL of Streptavidin C1 beads twice with 200uL of Streptavidin solution A, then 200uL of Streptavidin solution B. Discard the supernatants.
 2. Resuspend the beads in 25uL of 2X B&W buffer and add to the eluted SampleIp RNA. Incubate for 20-minutes at room temperature in a ThermoMixer with mixing at 1,000rpm.
 3. Wash the beads three times with 200uL of 1X B&W buffer and twice with 200uL of Streptavidin solution B. Discard the supernatants.
 4. Resuspend the beads in 10uL of water and incubate the samples at +90C for 5-minutes in a Thermocycler with mixing at 1,000rpm.
 5. Recover 8uL of the SampleIp RNA to a fresh tube and proceed to reverse transcription.
4. To perform reverse transcription for SampleIp and SampleInput:
1. In PCR tubes, add 2uL of 10uM Ribo_RTPrimer_v4, incubate at +70C for 5-minutes and then place on ice to cool.
 2. Add 10uL of the 2X RT master mix and pipette mix. Incubate in a Thermocycler at +55C for 30-minutes.
 3. Add 2.5uL of 10X Exo1 buffer followed by 2.5uL of Exonuclease I. Pipette mix and incubate at +37C for 15-minutes.
 4. Add 1uL of 0.5M EDTA to stop the reaction.
 5. Add 3uL of 1M NaOH and incubate at +70C in a Thermocycler for 12-minutes.
 6. Add 3uL of 1M HCl to neutralize the reaction.
 7. Magnetically separate 20uL of MyONE Silane beads per sample, and remove the supernatant.
 8. Rinse the beads once with 500uL of RLT buffer, and discard the supernatant.
 9. Reconstitute the beads in 370uL of RLT buffer and combine with the sample.
 10. Add 10uL of 5M NaCl to the sample and pipette mix.
 11. Add 400uL of ethanol and rotate the tubes to mix for 15-minutes at room temperature (+21C).
 12. Place tubes on magnetic rack and discard the supernatant.
 13. Rinse the beads three times with 800uL of 80% ethanol. Discard supernatants to the waste.
5. To ligate linkers to all samples:
1. Reconstitute the rinsed beads in 45uL of RNL buffer 2. Add 5uL of 3' linker mix and 2uL of T4 RNA ligase 1 (60 Units). Pipette mix and incubate at room temperature overnight with gentle agitation.
 2. Reconstitute the beads in 340uL of RLT buffer and combine with the sample.
 3. Add 10uL of 5M NaCl to the sample and pipette mix.
 4. Add 400uL of ethanol and rotate the tubes to mix for 15-minutes at room temperature (+21C).
 5. Place tubes on magnetic rack and discard the supernatant.

6. Rinse the beads three times with 800uL of 80% ethanol. Discard supernatants to the waste.
7. Add 25uL of water to the beads and incubate at room temperature for 5-minutes.
8. Recover the samples and store at -20C until PCR library construction.

5. Library amplification and size selection

5.1 Reagents and materials

- 1.5mL or 2.0mL Safe-Lock tubes (Thermo Scientific, CAT#05-402-25 or CAT#05-402-7)
- PCR microtubes (any vendor is fine)
- Benchtop centrifuge with holder for microcentrifuge tubes (multiple vendors)
- Centrifuge with buckets for microwell plates (multiple vendors)
- ThermoMixer with a 2mL block (Eppendorf)
- PCR Thermocycler (multiple vendors)
- RNA clean water (Thermo Scientific, CAT#10977023)
- SYBR Select Master Mix (Thermo Scientific, CAT#4472908)
- NEBNext Ultra II Q5 Master Mix (NEB, CAT#M0544L)
- MicroAmp Optical 384-Well Reaction Plate with Barcode (Thermo Scientific, CAT#4309849)
- MicroAmp Optical Adhesive Film (Thermo Scientific, CAT#4360954)
- Monarch PCR & DNA Cleanup Kit (5ug) (NEB, CAT#T1030S)
- SYBR Gold Nucleic Acid Gel Stain (10,000X Concentrate in DMSO) (Thermo Scientific, CAT#S11494)
- Novex Hi-Density TBE Sample Buffer (5X) (Thermo Scientific, CAT#LC6678)
- Novex TBE Running Buffer (5X) (Thermo Scientific, CAT#LC6675)
- Novex TBE Gels, 6%, 10 well (Thermo Scientific, CAT#EC6265BOX)
- pBR322 DNA-MspI Digest (NEB, CAT#N3032S)
- Sodium Acetate (3M), pH 5.5 (Thermo Scientific, CAT#AM9740)
- 5M NaCl (Thermo Scientific, CAT#AM9760G)
- 20% SDS (Thermo Scientific, CAT#BP1311)
- Disposable pellet pestles (VWR, CAT#KT749521-1500)
- Oligo clean and concentrator kit (Cedarlane, CAT#D4060)
- Razor blades (multiple vendors)
- 15mL tubes (VWR, CAT#89093-186)
- Absolute ethanol (wherever your institute purchases this from)
- Qubit dsDNA HS Assay Kit (Thermo Scientific, CAT#Q32851)
- Qubit assay tubes (Thermo Scientific, CAT#Q32856)
- DNA 1000 bioanalyzer kit (Agilent, CAT#5067-1504)

5.2 Solution recipes

- qPCR primer mix - 10uL of PCR_P5Primer_v4, 10uL of PCR_P7Primer_v4i1, 80uL water
- Gel ladder - mix 5uL of pBR322 DNA-MspI Digest, 5uL of Novex Hi-Density TBE Sample Buffer, and 15uL of water
- Gel stain - 5uL of SYBR Gold in 50mL of 1X TBE
- Gel elution buffer (make 10mL)
 - 300mM NaCl (600uL of 5M stock)
 - 60mM sodium acetate pH 5.5 (200uL of 3M stock)
 - 0.2% SDS (100uL of 20% SDS)
 - water to 10mL (9.1mL)

5.3 Oligos

- PCR_P5Primer_v4 (4nmol scale, Ultramer)
 - /5Phos/AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATC
 - reconstitute to 100uM and store at -20C.
- PCR_P7Primer_v4i1 (4nmol scale, Ultramer)
 - /5Phos/CAAGCAGAAGACGGCATACGAGATCGTGATGTGACTGGAGTTCAGACGTGTGCTCTA
 - reconstitute to 100uM and store at -20C.
- PCR_P7Primer_v4i2 (4nmol scale, Ultramer)
 - /5Phos/CAAGCAGAAGACGGCATACGAGATACATCGGTGACTGGAGTTCAGACGTGTGCTCTA
 - reconstitute to 100uM and store at -20C.
- PCR_P7Primer_v4i3 (4nmol scale, Ultramer)
 - /5Phos/CAAGCAGAAGACGGCATACGAGATGCCTAAGTGACTGGAGTTCAGACGTGTGCTCTA
 - reconstitute to 100uM and store at -20C.
- PCR_P7Primer_v4i4 (4nmol scale, Ultramer)
 - /5Phos/CAAGCAGAAGACGGCATACGAGATTGGTCAGTGACTGGAGTTCAGACGTGTGCTCTA
 - reconstitute to 100uM and store at -20C.
- PCR_P7Primer_v4i5 (4nmol scale, Ultramer)
 - /5Phos/CAAGCAGAAGACGGCATACGAGATCACTGTGTGACTGGAGTTCAGACGTGTGCTCTA
 - reconstitute to 100uM and store at -20C.
- PCR_P7Primer_v4i6 (4nmol scale, Ultramer)
 - /5Phos/CAAGCAGAAGACGGCATACGAGATATTGGCGTGACTGGAGTTCAGACGTGTGCTCTA
 - reconstitute to 100uM and store at -20C.
- PCR_P7Primer_v4i7 (4nmol scale, Ultramer)
 - /5Phos/CAAGCAGAAGACGGCATACGAGATGATCTGGTGACTGGAGTTCAGACGTGTGCTCTA
 - reconstitute to 100uM and store at -20C.
- PCR_P7Primer_v4i8 (4nmol scale, Ultramer)
 - /5Phos/CAAGCAGAAGACGGCATACGAGATTCAAGTGTGACTGGAGTTCAGACGTGTGCTCTA
 - reconstitute to 100uM and store at -20C.
- PCR_P7Primer_v4i9 (4nmol scale, Ultramer)
 - /5Phos/CAAGCAGAAGACGGCATACGAGATCTGATCGTGACTGGAGTTCAGACGTGTGCTCTA
 - reconstitute to 100uM and store at -20C.
- PCR_P7Primer_v4i10 (4nmol scale, Ultramer)

- /5Phos/CAAGCAGAAGACGGCATAACGAGATAAGCTAGTGACTGGAGTTCAGACGTGTGCTC
- reconstitute to 100uM and store at -20C.
- PCR_P7Primer_v4i11 (4nmol scale, Ultramer)
 - /5Phos/CAAGCAGAAGACGGCATAACGAGATGTAGCCGTGACTGGAGTTCAGACGTGTGCTC
 - reconstitute to 100uM and store at -20C.
- PCR_P7Primer_v4i12 (4nmol scale, Ultramer)
 - /5Phos/CAAGCAGAAGACGGCATAACGAGATTACAAGGTGACTGGAGTTCAGACGTGTGCTC
 - reconstitute to 100uM and store at -20C.

5.4 Protocol

The starting point for this stage of the protocol is the 25uL purified cDNA samples from above.

1. First, determine the amplification levels of the different libraries:
 1. Combine 1uL of cDNA sample, 1uL of qPCR primer mix, 8uL of water, and 10uL of SYBR Select master mix and pipette mix.
 2. Add sample to 384-well qPCR plate and seal with optical adhesive film. Centrifuge at 250g for 5-minutes.
 3. Analyze on an RT-qPCR instrument with using the cycle settings:
 1. +95C for 10-seconds
 2. 40-cycles of:
 3. +95C for 15-seconds
 4. +60C for 30-seconds
 5. +72C for 30-seconds
 6. Melt curve analysis
 7. +4C forever
 4. Note the Ct values for the samples (SampleIp optimal cycles = Ct - 4, SampleInput optimal cycles = Ct - 1).
2. Perform a test amplification to see that the insert size is correct:
 1. In a PCR tube, combine 1uL of sample, 1uL of 10uM qPCR primer mix, 8uL of water, and 10uL of Q5 PCR master mix and pipette mix.
 2. Amplify in a Thermocycler with the cycle settings:
 1. +98C for 30-seconds
 2. X-cycles of:
 3. +98C for 10-seconds
 4. +65C for 75-seconds
 5. +65C for 5-minutes
 6. +4C forever
 3. Add 5uL of 5X Novex Hi-Density TBE Sample Buffer to the PCR reaction and run on a TBE gel.
 4. Stain in Gel stain and visualize the result.
3. Perform library amplification (SampleIp optimal cycles = Ct - 6, SampleInput optimal cycles = Ct - 2).:
 1. In a PCR tube, combine 12.5uL of sample, 2.5uL of 10uM primer mix (with your barcode of interest), and 12.5uL of Q5 PCR master mix

- and pipette mix.
2. Amplify in a Thermocycler with the cycle settings:
 1. +98C for 30-seconds
 2. X-cycles of:
 3. +98C for 10-seconds
 4. +65C for 75-seconds
 5. +65C for 5-minutes
 6. +4C forever
 3. Purify the PCR product with a Monarch PCR purification kit and elute in 20uL of water.
 4. Add 5uL of 5X Novex Hi-Density TBE Sample Buffer to the PCR reaction and run on a TBE gel (120V). Use 5uL of the Gel ladder per lane as a marker.
 5. Stain Gel stain and visualize the result.
 6. Excise the region between 175 to 350bp and place into a fresh 1.5mL tube.
 7. Using a pellet pestle, crush the gel piece in the tube.
 8. Add 400uL of Gel elution buffer and freeze the tube on dry ice. Rotate the tube at room temperature for at least 2-hours.
 9. Centrifuge the tube at 12,000g for 2-minutes and transfer the supernatant to a 15mL tube.
 10. From the Oligo clean and concentrator kit, add 800uL of Oligo binding buffer to the sample (some precipitation may occur here, that is ok).
 11. Add 1.6mL of ethanol and vortex mix.
 12. Transfer 750uL of sample to an Oligo clean and concentrator column and centrifuge at 10,000g for 30-seconds. Discard the flow-through.
 13. Repeat the previous centrifuge step until all sample material has been processed through the column.
 14. Add 750uL of DNA wash buffer to the column and centrifuge for 30-seconds at 10,000g. Discard the flow-through.
 15. Centrifuge the column for 1-minute at 12,000g.
 16. Transfer the column to a fresh 1.5mL collection tube, add 15uL of water to the filter, and incubate for 1-minute at room temperature.
 17. Centrifuge for 1-minute at 10,000g to recover the sample material.
 4. Measure the concentration of 1uL of the sample using the Qubit HS DNA kit.
 5. Analyze the library with a Bioanalyzer DNA kit.
 6. If multiplexing, combine the library based on the Qubit concentrations.
 7. Submit for sequencing (libraries are compatible with single-end + index read on NextSeq, or paired end on all Illumina systems).

Notes

Note 1 - This setting seems to be different in every paper you read about CLIP. You may want to test the efficiency here yourself, but in general I find that most publications fall in a range of 150 - 400mJ/cm². Some studies have reported less

noise in the sequencing data at lower crosslinking energies, so I tend towards this side of that range. There is some nice additional discussion of crosslinking parameters in the paper PMID: 24184352. When crosslinking, the lid should be off your tissue culture dish. For the ice block, I like to use the lid off of a 15cm dish where I have added a small amount of water to it and frozen it. In the crosslinker, I will add a small amount of water on top so the cell culture dish doesn't have to sit directly on the ice (this is adapted from the irCLIP2 protocol).

Note 2 - This antibody concentration is dependent on your own antibody. You should validate the specificity and efficiency of the antibody you plan to use. An easy way to do this is to do some test pull-downs and analyze the products via SDS-PAGE, Western, or mass spectrometry.

Note 3 - The addition of BSA here serves to reduce non-specific binding on the beads and antibody. This step can be skipped if you are not concerned about this.

Note 4 - As written, this protocol will also isolate RNA from nuclei. If you wish to only look at cytoplasmic RNA, simply exclude SDS from the Lysis buffer and skip this mechanical disruption step. Reduce the speed of the debris centrifugation step (the next step) from 12,000g to 5,000g. The pellet will contain nuclei. The mechanical disruption step can be excluded entirely if desired, as most RNA should leak out of the nucleus in the presence of 0.1% SDS.

Note 5 - To recover more of the liquid here, you can use a small needle that is heated in a bunsen burner flame to poke a hole in the bottom of the bead-containing tube. Insert the tube with the hole into a standard 2mL tube and spin at 50g to draw the liquid through the hole into the 2mL tube. Alternatively, you can just sonicate here.

Note 6 - The RNase 1 concentration here is very important. Generally, you want the amount of RNase you are adding to each of your samples to be equal based on how much RNA they have. I find it is best to normalize the lysate concentrations and go from there. I have seen it suggested to use Nanodrop assays here, but this is risky because it is so inaccurate and prone to interference, especially when working with a whole-cell lysate. A BCA assay would theoretically be ok, but it doesn't tell us anything about the RNA amount in the sample, which could change in certain conditions independent of protein levels. Instead of relying on a 'dilution ratio' like so many studies, I like to provide an exact number of Units of RNase to add based on how much RNA we have in each sample.

Note 7 - Important to remember that RNase 1 still has activity at +4C, so be mindful of how long your incubation time is here and do not vary it between samples that are to be directly compared.

Reagent list

I know it can be a bit hard to parse all of the reagents above, especially since many are repeated in the different sections. To make this easier, this is a list of the required reagents in an easier to parse format. It excludes things like tubes and equipment pieces. Prices are taken directly from vendor websites and are current as of May 2020, in Canadian dollars.

ItemNumber	Description	Vendor	Quantity	Price
AM2238	TURBO Dnase, 2U/uL, 1000 Units	Thermo Scientific	1	\$152.00
N2511	RNaseIn, 40U/uL, 2500 Units	Promega	1	\$139.00
15567027	UltraPure 1M TrisCl, pH 7.5, 1L	Thermo Scientific	1	\$75.75
AM9760G	NaCl, 5M, Rnase free, 500mL	Thermo Scientific	1	\$108.00
15575020	UltraPure 0.5M EDTA, pH 8, 4x100mL	Thermo Scientific	1	\$109.00
NP40	IGEPAL CA-630 solution, 100mL	Sigma	1	\$115.00
D6750-100G	Sodium deoxycholate, 100G	Sigma	1	\$150.00
AM9530G	MgCl ₂ , 1M, 100mL	Thermo Scientific	1	\$50.75
11836170001	cOmplete protease inhibitor tablets, EDTA free, 25 tablets	Sigma	1	\$202.00
Q32856	Qubit assay tubes, 100 tubes	Thermo Scientific	1	\$120.00

ItemNumber	Description	Vendor	Quantity	Price
Q32852	Qubit RNA HS assay kit, 100 assays	Thermo Scientific	1	\$144.00
AM2294	RNaseI, cloned, 100U/uL, 10,000 Units	Thermo Scientific	1	\$188.00
10977023	UltraPure Dnase/Rnase free water, 10x500mL	Thermo Scientific	1	\$232.00
116960050	Lysing matrix Y tubes, 50x2mL	Cedarlane	1	\$299.00
NP0007	NuPAGE LDS Sample Buffer (4X), 10mL	Thermo Scientific	1	\$34.87
15503022	UltraPure Sucrose, 5kg	Thermo Scientific	1	\$249.00
BP1311	20% SDS, Nuclease free, 1L	Thermo Scientific	1	\$148.56
88802	Protein A/G magnetic beads, 1mL	Thermo Scientific	1	\$392.00
A2153-10G	Bovine serum albumin, 10G	Sigma	1	\$122.00
1610611	Dithiothreitol, 5g	Bio-Rad	1	\$121.00
V3011	PEG-8000, 500g	Promega	1	\$87.00
M0525S	QuickCIP, 1000 Units	NEB	1	\$131.00
M0201L	T4 PNK, 2500 Units	NEB	1	\$297.00
M0373L	T4 RNA ligase 2, truncated KQ, 10,000 Units	NEB	1	\$362.00
P0756S	10mM ATP, 1mL	NEB	1	\$46.80

ItemNumber	Description	Vendor	Quantity	Price
M0437M	T4 RNA ligase 1, high concentra- tion, 5000 Units	NEB	1	\$335.00
26616	PageRuler prestained protein ladder, 10 to 180kDa, 2x250uL	Thermo Scientific	1	\$204.00
NP0335BOX	NuPAGE 4 to 12% Bis- Tris, 1.5 mm, Mini-Scientific tein Gel, 10- well	Thermo Scientific	1	\$265.00
NP0001	NuPAGE MOPSDS Run- ning Buffer (20X), Scientific 500mL	Thermo Scientific	1	\$114.00
NP00061	NuPAGE Trans-Thermo fer Buffer (20X), Scientific 1L	Thermo Scientific	1	\$349.00
GE10600008	Amersham Pro- Sigma tran Pre- mium West- ern blot- ting mem- branes, nitrocellulose		1	\$376.00
88518	PVDF Trans- Thermo fer Mem- Scientific brane, 0.45um		1	\$618.00
34860-1L-R	Methanol, Sigma HPLC grade, 1L		1	\$67.90
18-8816-31	TrueBlot Anti- Cedarlane Rabbit IgG HRP		1	\$233.00
14190144	DPBS, no Thermo calcium, no Scientific magnesium, 500mL		1	\$40.51
P9416- 100ML	Tween 20, Sigma 100mL		1	\$44.60
93362-250G	Trizma base, Sigma 250G		1	\$184.00

ItemNumber	Description	Vendor	Quantity	Price
71376-1KG	Sodium chloride, 1KG	Sigma	1	\$140.00
H1758-100ML	Hydrochloric acid, 100mL	Sigma	1	\$84.20
32109	Pierce ECL Western Blotting Substrate, 50mL	Thermo Scientific	1	\$101.00
34089	CL-XPose Film, 100 sheets	Thermo Scientific	1	\$338.00
AM2546	Proteinase K Solution (20 mg/mL), 1.25mL	Thermo Scientific	1	\$91.50
U5378-500G	Urea, 500G	Sigma	1	\$89.10
AM9720	Acid-Phenol:Chloroform:isoamyl alcohol, 100mL	Thermo Scientific	1	\$116.00
P3803-100ML	Phe-nol:Chloroform:Isoamyl Alcohol 25:24:1, Saturated with 10mM Tris, pH 8.0, 1mM EDTA	Sigma	1	\$129.00
CA10847-802	Phaselock gel (heavy), 200 tubes	VWR	1	\$533.75
R1015	RNA Clean and concentrator-5, 50 preps	Cedarlane	1	\$220.00
37002D	Dynabeads MyOne Silane, 5mL	Thermo Scientific	1	\$495.00
79216	RLT buffer, 220mL	Qiagen	1	\$167.00
65001	Dynabeads MyOne Streptavidin C1, 2mL	Thermo Scientific	1	\$724.00

ItemNumber	Description	Vendor	Quantity	Price
M0293S	Exonuclease I, 2000 Units (E. coli), 3000 Units	NEB	1	\$90.90
S2770-100ML	1N NaOH, 100mL	Sigma	1	\$24.70
H9892-100ML	1N HCl, 100mL	Sigma	1	\$32.80
18090050	SuperScript IV Reverse Transcriptase, 10,000 Units	Thermo Scientific	1	\$508.00
N0447S	Deoxynucleotide (dNTP) Solution Mix, 10mM	NEB	1	\$81.00
4472908	SYBR Select master mix, 5mL	Thermo Scientific	1	\$500.00
M0544L	NEBNext UltraII Q5 master mix, 50rxns	NEB	1	\$131.00
4309849	MicroAmp Optical 384-Well Reaction Plate with Barcode, 50 plates	Thermo Scientific	1	\$485.00
4360954	MicroAmp Optical Adhesive Film, 25 sheets	Thermo Scientific	1	\$115.00
T1030S	Monarch PCR & DNA cleanup kit, 50 preps	NEB	1	\$115.20
S11494	SYBR Gold nucleic acid gel stain, 500uL	Thermo Scientific	1	\$224.00
LC6678	Novex Hi-Density TBE sample buffer, 5X, 10mL	Thermo Scientific	1	\$147.00

ItemNumber	Description	Vendor	Quantity	Price
LC6675	Novex TBE Running buffer, 5X, 1L	Thermo Scientific	1	\$62.00
EC6265BOX	Novex TBE Gels, 6%, 10-well, box of 10	Thermo Scientific	1	\$252.00
N3032S	pBR322 DNA-MspI Digest, 50uL	NEB	1	\$96.30
AM9740	Sodium Acetate (3M), pH 5.2, 100mL	Thermo Scientific	1	\$50.75
KT749521-1500	Disposable pellet pestles, 100 pack	VWR	1	\$243.59
D4060	Oligo clean and concentrator kit, 50 preps	Cedarlane	1	\$133.00
Q32851	Qubit ds-DNA HS Assay Kit, 100 assays	Thermo Scientific	1	\$143.00
5067-1504	Agilent DNA 1000 kit, 25 chips	Agilent	1	\$1,267.00